

Cloning of DNA Fragments

Created: December 27, 1994
Last Revision: October 24, 2012

A. Ligation

1. Purify vector (e.g., pUC19, pET11d) by commercial plasmid prep (e.g., Quiagen QUIAprep Spin Plasmid Kit) and run on an agarose gel to determine purity and concentration. Additionally, determine OD₂₆₀ and OD₂₈₀ of a 1:100 dilution in 250 µl cuvetts in UV spectrophotometer.
2. Digest 1 to 2 µg (or as much as possible) of plasmid in a final volume of 40 µl with the appropriate restriction enzymes. Purify plasmid DNA by phenol and chloroform extraction followed by ethanol precipitation with 1/10 volume NaCl. Dissolve pellet in 9 µl dest. water and check restriction digest and DNA concentration by running 2 µl on an agarose gel.
3. Digest 1 to 2 µg (or as much as possible) of the target DNA with the appropriate enzymes. Purify DNA after restriction digest. For PCR products, use Quiagen Kit and dry DNA afterwards in vacuum centrifuge (SpeedVac). Otherwise, purify by phenol and chloroform extraction followed by ethanol precipitation with 1/10 volume NaCl. Dissolve pellet in 9 µl dest. water and check restriction digest and DNA concentration by running 2 µl on an agarose gel.
4. Set up ligation in a 20 µl volume with a vector:insert ratio of cohesive ends between 1:1 and 3:1. Set up a negative control ligation with vector and no insert.

digested and purified vector DNA:	100 ng
digested and purified insert DNA:	1 to 3 times as many cohesive ends as vector
10X ligase buffer:	2 µl
T4 ligase:	approx. 0.1 Weiss Unit
Dest. water:	up to 100 µl

Incubate @ 15°C overnight.

B. Transformation

1. Mix about 10 µl of the ligation and 200 µl competent cells (prepared as described in separate protocol), keep on ice for 45 min.
2. Transfer cells into tube preheated @ 42°C and incubate @ 42°C for **exactly** 90 sec.; chill in ice bath for 1 - 2 min.

3. Add 800 μ l LB to cells and incubate @ 37°C (on shaker) for 45 min.
4. Plate cells on LB/X-gal/IPTG (for blue white screening) or on LB plates with appropriate antibiotic (e.g., 100 μ g/ml ampicillin) and incubate @ 37°C for 12 to 24 h. Plate 10 μ l, 50 μ l and 200 μ l of the transformation. Keep remainder of transformation at 4°C.

References:

Sambrook J., E. F. Fritsch, and T. Maniatis (1989) Molecular cloning: a laboratory manual, 2nd ed. Cold Spring Harbor Laboratory, Cold Spring Harbor, N.Y.