



FOOD SAFETY LAB / MILK QUALITY IMPROVEMENT PROGRAM  
*Standard Operating Procedure*

Title: **ITS Yeast and Mold PCR Protocol**

SOP #:

Revision: **00**

Revision Date: **n/a**

Effective Date: **2019-11-14**

Authors: **Jordan Skeens, Nicole Martin, Rachel Miller**

Approved by:

***ITS Yeast and Mold PCR Protocol***

**FILE NAME: ITS Yeast and Mold PCR Protocol\_10112019.doc**



## TABLE OF CONTENTS

1. INTRODUCTION
  - Purpose
  - Scope
  - Definitions
  - Safety
2. MATERIALS
3. PROCEDURE
4. TROUBLESHOOTING
5. REFERENCES



## **SECTION 1      INTRODUCTION**

### **1.1      Purpose**

This protocol describes the standardized method to amplify a ~700 bp portion of the ITS gene for yeast and mold by PCR for sequencing and subsequent organism identification.

This SOP is derived from an earlier protocol titled “Yeast and Mold PCR SOP”. The original protocol can be found in the archived section of the Food Safety Wiki.

### **1.2      Scope**

This SOP applies to the Food Safety Lab, including the Milk Quality Improvement Program. The protocols may also be used by laboratory members from other locations.

### **1.3      Definitions**

**BP:** base pair

**dNTPs:** deoxyribonucleotide monomers, which are the building blocks of DNA. There are four dNTPs, cytosine, adenine, guanine and thymine

**GoTaq Green:** A PCR kit containing thermostable Taq polymerase, MgCl<sub>2</sub>, and a PCR buffer which contains a green loading dye. This loading dye does not interfere with the PCR and saves the step of loading dye incorporation prior to running gel electrophoresis.

**PCR:** polymerase chain reaction, used to amplify a specific region within a DNA sequence

### **1.4      Safety**

Wear gloves for the entire procedure.



## SECTION 2 MATERIALS

### Single Tube Reactions

- Sterile water - Room 358B and Room 350A
- Sterile 0.2 mL tubes - Room 358B and Room 350A
- ITS 4 (forward) Primer, 10 mM (5' - TCC TCC GCT TAT TGA TAT GC -3') -Room 352C, chest freezer
- ITS 5 (reverse) Primer 10mM (5' - GGA AGT AAA AGT CGT AAC AAG G - 3') - Room 352C, chest freezer
- GoTaq Flexi Polymerase - Room 352, “Dumb” freezer
- GoTaq Green 5x PCR Buffer - Room 352, “Dumb” freezer
- dNTP solution consisting of 10 mM each dATP, dGTP, dCTP, dTTP dNTP stocks in Room 352, “Dumber” freezer, in covered cryoblock
- MgCl<sub>2</sub>, 25 mM - Room 352, “Dumb” freezer
- Vortex Mixer - Room 358B
- Miniature Centrifuge - Room 358B
- Thermocycler - Room 356
- Micropipette and sterile filter tips
- Crushed Ice

### High Throughput Reactions (96-well plate)

- Sterile reservoir basin - Room 358B
- Sterile 96-well plates - Room 358B
- Multi-channel pipette
- Sterile water - Room 358B and Room 350A
- ITS 4 (forward) Primer, 10 mM (5' - TCC TCC GCT TAT TGA TAT GC -3') -Room 352C, chest freezer
- ITS 5 (reverse) Primer 10mM (5' - GGA AGT AAA AGT CGT AAC AAG G - 3') - Room 352C, chest freezer
- GoTaq Flexi Polymerase - Room 352, “Dumb” freezer
- GoTaq Green 5x PCR Buffer - Room 352, “Dumb” freezer
- dNTP solution consisting of 10 mM each dATP, dGTP, dCTP, dTTP. Stocks in Room 352, “Dumber” freezer, in covered cryoblock
- MgCl<sub>2</sub>, 25 mM - Room 352, “Dumb” freezer
- Vortex Mixer - Room 358B
- Miniature Centrifuge - Room 358B
- Thermocycler - Room 356
- Adhesive aluminum foil lids - Room 358B
- Micropipette and sterile filter tips
- Crushed Ice



## SECTION 3 PROCEDURES

### 3.1 Lysate Preparation

Refer to the “Preparing DNA Lysates” SOP on the Food Safety Lab wiki.

### 3.2 PCR Amplification

- 3.2.1 Working in Room 358B, prepare a master mix of the components listed in Table 1 by pipetting each component into a sterile 1.5 mL microcentrifuge tube.
- 3.2.1.1 Pipette up and down to mix after each addition and to ensure that no reagent clings to the pipette tip.
- 3.2.1.2 Store all master mix components and the prepared master mix on crushed ice.
- 3.2.1.3 Vortex and centrifuge down the contents of all reagent containers before opening.
- 3.2.1.4 Calculate master mix volumes required to include positive and negative controls and 1-2 additional reactions to compensate for pipetting errors.

**Table 1: ITS PCR Master Mix**

<u>Reagent</u>	<u>Volume per reaction (µL)</u>
dH2O	32.5 µl
5X Green GoTaq Flexi Buffer	10 µl
dNTPs 1mM	1 µl
MgCl <sub>2</sub> 25mM	2 µl
ITS 4 Primer 10mM	1 µl
ITS 5 Primer 10mM	1 µl
GoTaq Flexi DNA polymerase	0.5 µl
<hr/>	
<b>Total</b>	<b>48.0 µL</b>

- 3.2.2 Briefly vortex the prepared master mix. Spin briefly in the mini centrifuge to concentrate all contents at the bottom of the tube.



- 3.2.3 Dispense 48  $\mu$ L of prepared master mix into one well per lysate in a 96-well PCR plate or into individual 0.2 mL PCR tubes. Cover the plate loosely or cover the tubes and hold on ice.
- 3.2.4 Remove the plate or tubes containing the master mix from room 358B to a lab bench.
- 3.2.5 Briefly spin down the previously-prepared lysates in a centrifuge, then add 2  $\mu$ L of lysate to its corresponding PCR well or tube. Pipette up and down to mix.
- 3.2.6 Briefly centrifuge the PCR plate or tubes to concentrate reagents at the bottom of the tube.
- 3.2.7 Set up the thermocycler using the parameters described in 3.2.7.1.

**3.2.7.1 Thermocycling Conditions for ITS PCR**

Time (minutes:seconds)	Temperature ( $^{\circ}$ C)	Number of Cycles
5:00	95	1
1:00	95	35
1:00	56	
1:00	72	
10:00	72	1
$\infty$	4	1

- 3.2.8 Promptly remove the reactions from the thermocycler when complete. PCR product can be stored at 4 $^{\circ}$ C for up to 1 week or at -20 $^{\circ}$ C for up to 6 months.

**3.3 Gel Electrophoresis Confirmation of Amplification**

Refer to the “Gel Electrophoresis of PCR Product” SOP on the Food Safety Lab wiki.

**3.4 ExoSAP Purification and Sequencing Submission**

Refer to the “PCR Product Purification and Sanger Sequencing Submission” SOP on the Food Safety Lab wiki.

**3.5 Sequence Analysis**

Refer to the “Automated Sequence Editing” SOP on the Food Safety Lab wiki.



## SECTION 4 TROUBLESHOOTING

### 4.1 Issues with MgCl<sub>2</sub>

- 4.1.1 MgCl<sub>2</sub> forms a concentration gradient when frozen and needs to be vortexed prior to use.
- 4.1.2 Every PCR reaction has an optimal MgCl<sub>2</sub> concentration range, usually between 1-4 mM. Since Mg<sup>+2</sup> forms complexes with dNTPs and acts as a cofactor for polymerase you may need to try several conditions to optimize your concentration.

### 4.2 Too Much Enzyme

- 4.2.1 Excess enzyme in your PCR can lead to smearing of PCR products.

### 4.3 Wrong Primer Concentration

- 4.3.1 If you have too little primer you won't see any product.
- 4.3.2 Too much primer may result in primer dimerization, which may reduce amplification.

### 4.4 Wrong PCR Program

- 4.4.1 Check your program while it's cycling to make sure it is the right program.

### 4.5 Excess or Insufficient Template

- 4.5.1 Too much template can inhibit PCR by binding all of the primers. If the lysate appears very turbid, dilute the template 1:10 or 1:100 and then use 1 µL of the diluted template for the PCR.
- 4.5.2 Too little template and amplification may be undetectable. Make sure you correctly prepared your dirty lysate and that you took a sufficiently large sample from your culture.

### 4.6 Too Much dNTP or Degraded dNTP

- 4.6.1 Excess dNTP inhibits PCR due to Mg<sup>2+</sup> depletion. Between 40-200 µM is optimal. dNTPs are sensitive to repeated freeze-thaw cycles. Make small aliquots when you get a fresh batch and turn over your stock frequently.
- 4.6.2 Note that the dNTP solution concentration suggested by the GoTaq Flexi manufacturer Promega suggests having each dNTP at 10mM, while this protocol uses 1mM due to less Mg<sup>2+</sup> being used.

### 4.7 Suggested Positive control: FSL E2-0320 (*Penicillium decumbens*), stored in freezer 3, tower 25. Amplifies with both yeast and mold primer sets.



4.8 The NL gene can be amplified for PCR troubleshooting if the ITS gene fails to amplify. Use the same reagents and volumes as listed in Materials but use the NL primer set and thermocycle settings:

- NL 1 (forward) Primer, 10 mM (5' - ATA TCA ATA AGC GGA GGA AAA G -3') - Room 352C, chest freezer
- NL 4 (reverse) Primer 10mM (5' - GGT CCG TGT TTC AAG ACG G - 3') - Room 352C, chest freezer

**Thermocycling Conditions for NL PCR**

<b>Time (minutes:seconds)</b>	<b>Temperature (°C)</b>	<b>Number of Cycles</b>
5:00	94	1
1:00	94	36
1:00	52	
2:00	72	
10:00	72	1
∞	4	1





## **SECTION 5**

## **REFERENCES**

Yeast and Mold DNA Extraction SOP

Primer Dilution SOP

dNTPs PCR Nucleotide Mix SOP