

## Sybr Green Real Time qRT-PCR

**Purpose:** To assay transcript levels using SYBR® Green.

### Reagents:

- SYBR® Green PCR Master Mix (Applied Biosystems®) Catalog Number 4309155 (Size 1 x 5 mL)
- TaqMan® Reverse Transcription Reagents (Invitrogen™) Catalog Number N8080234 (Size 200 reactions)
- MicroAmp® Optical 96-Well Reaction Plate (Applied Biosystems®) Catalog Number N8010560 (Size 10 plates)
- MicroAmp® Optical Adhesive Film (Applied Biosystems®) Catalog Number 4360954 (Size 25 covers).

### Protocol:

#### Step 1. Reverse Transcription (From TaqMan® Universal PCR Master Mix)

A 100- $\mu$ L RT reaction efficiently converts a maximum of 2  $\mu$ g total RNA to cDNA. Perform multiple RT reactions in multiple wells if you are using more than 2  $\mu$ g total RNA. For a 10- $\mu$ L RT reaction you can use up to 500ng of RNA in  $\leq$  3.85 $\mu$ l of water.

#### RT Reaction Mix

Component	Volume ( $\mu$ L) Per Sample	Final Concentration
RNase-free water	See below <sup>a</sup>	—
10X TaqMan RT Buffer	1.0	1X
25 mM Magnesium Chloride	2.2	5.5 mM
deoxyNTPs Mixture	2.0	500 $\mu$ M per dNTP
Random Hexamers <sup>b</sup>	0.5	2.5 $\mu$ M
RNase Inhibitor	0.2	0.4 U/ $\mu$ L
MultiScribe Reverse Transcriptase (50 U/ $\mu$ L)	0.25	1.25 U/ $\mu$ L
<b>Total</b>	<b>6.15<sup>c</sup></b>	<b>—</b>

a. The volume of RNase-free water ( $\mu$ L) is 3.85–RNA sample volume in a 10- $\mu$ L reaction.

b. Random hexamers, oligo d(T)<sub>16</sub>, or sequence-specific reverse primers can be used for primers of cDNA synthesis.

c. If changing the reaction volume, make sure the final proportions are consistent with the recommended values above.

**Note** The RT volume can vary from 10  $\mu$ L to 100  $\mu$ L. Increasing the RT volume will reduce the total number of reactions.

**Thermal Cycling** To conduct RT thermal cycling:

Step	Action																				
1	Load the reactions into a thermal cycler.																				
2	Program the thermal cycler or Sequence Detection System with the following conditions: <table border="1" data-bbox="532 409 1205 661"> <thead> <tr> <th>Step</th> <th>Incubation<sup>a</sup></th> <th>RT</th> <th>Reverse Transcriptase Inactivation</th> </tr> </thead> <tbody> <tr> <td></td> <td>HOLD</td> <td>HOLD</td> <td>HOLD</td> </tr> <tr> <td>Temp</td> <td>25 °C</td> <td>48 °C</td> <td>95 °C</td> </tr> <tr> <td>Time</td> <td>10 min</td> <td>30 min</td> <td>5 min</td> </tr> <tr> <td>Volume</td> <td colspan="3">10 µL</td> </tr> </tbody> </table> <p>a. If using random hexamers or oligo d(T)<sub>16</sub> primers for first-strand cDNA synthesis, a primer incubation step (25 °C for 10 min) is necessary to maximize primer-RNA template binding.</p>	Step	Incubation <sup>a</sup>	RT	Reverse Transcriptase Inactivation		HOLD	HOLD	HOLD	Temp	25 °C	48 °C	95 °C	Time	10 min	30 min	5 min	Volume	10 µL		
Step	Incubation <sup>a</sup>	RT	Reverse Transcriptase Inactivation																		
	HOLD	HOLD	HOLD																		
Temp	25 °C	48 °C	95 °C																		
Time	10 min	30 min	5 min																		
Volume	10 µL																				
3	Begin RT.  <b>IMPORTANT</b> After thermal cycling, store all cDNA samples at -15 to -25 °C.																				

### Step 2. qRT-PCR using sybr green.

- Dilute your cDNA 1/10 in water.
- Prepare master mix for each primer set. Final volume per well= 25 µl.

Reagent	Volume
2X SYBR	12.5 µl
Water	10.5 µl
10µM Primer 1	0.5 µl
10µM Primer 2	0.5 µl
Template (1-5ng)	1 µl

### Step 3. Assay and Analysis.

#### A. Load Optical 96 Well Reaction Plate

- Pipet 1 µl from each template into appropriate well.
- Pipet 24 µl from master mix into appropriate well.
- Load 25 µl from leftover Master Mix tubes as No-Template-Controls (NTC). Load at least one NTC for each primer pair.
- Cover the plate with optical adhesive cover film. Handle only by the edges and do not touch the surface of the film

- Spin the plate for 2 min. in the benchtop swing bucket centrifuge (1000 rpm for one minute).
- Run standard Real Time PCR in the QuantStudio 6 system.