

Protocol for Passaging HIEC-6 cells

Materials

Pasteur Pipettes

Flasks (T25 or T75)

Water Bath at 37C

Stripettes (2 mL 5 mL, 10 mL)

PBS

0.25% Trypsin EDTA (ThermoFisher Scientific catalog #: 25200072)

Opti-MEM media (can be stored up to 2 weeks at 4C):

- 10 ng/mL EGF (ThermoFisher Scientific catalog #: PHG0315)
- Opti-MEM supplemented with HEPES and Glutamax (Gibco catalog #: 51985034)
- 10% FBS
- *I have cultured the cells with both 4% FBS (as recommended on ATCC website) and 10% FBS in parallel, but I did not notice a difference.*

37C Incubator with 5% CO₂

1) Preparation

- Turn on 37C water bath. Let temperature equilibrate to 37°C (approx. 15 mins).
- Place the PBS and HIEC-6 media in the 37°C water bath. Allow temperature to equilibrate to 37°C (approx. 15 mins).
- Allow Trypsin solution to warm to room temperature (15 mins).
- Turn on UV light in biosafety cabinet (approx. 15 mins)
- After 15 mins, turn off UV light in the biosafety cabinet. Turn on the fluorescent light and blower. Open the biosafety cabinet and spray down surface with ethanol (liberally!).
- Transfer all materials to the hood (stripettes, PBS, Trypsin, Opti-MEM, new flasks)

2) Trypsinize cells

- Remove flask of HIEC-6 cells from the 37°C incubator. Check with inverted microscope to make sure that flask is between 90-100% confluent.
- In biosafety cabinet, remove media with a Pasteur pipette and vacuum pump.
- Rinse cells by adding 2-5 mL of pre-warmed PBS to the flask. Remove with a Pasteur pipette.
- Add 0.5 – 1mL (0.5 mL for T25 flask, 1 mL for T75 flask) of room temperature trypsin-EDTA. Return flask to the 37C incubator for 5 minutes. Check for evidence of detachment (cells will appear as “flakes”). *Note: HIEC-6 cells require longer to detach (compared to HeLa cells). I have not noticed any problems incubating them with Trypsin-EDTA for >5 mins.*
- Stop trypsinization by adding Opti-MEM (4.5 mL for T25 and 9 mL for T75). Pipette cells up and down 8-10 times to break up clumps.

3) Seed a new flask.

- For HeLa cells, we will use the following split ratios. Add the appropriate amount of cell suspension and media, and transfer to a new flask.

<u>Split Ratio</u>	<u>mL cells /mL media for T25</u>	<u>mL cells /mL media for T75</u>
1:4	1 mL + 4 mL	2 mL + 8 mL
1:5	0.8 mL + 4.2 mL	1.6 mL + 8.4 mL
1:6	0.7 mL + 4.3 mL	1.4 mL + 8.6 mL
1:8	0.6 mL + 4.4 mL	1.2 mL + 8.8 mL
1:10	0.5 mL + 4.5 mL	1.0 mL + 9 mL

- Label flask! (cell line, date, passage #, initials) and return to the 37°C incubator. *These cells do not do well if the split ratio is too high. 1:4 or 1:5 is optimum.*

For Generating Freezer Stocks

- Retrieve pre-chilled isopropanol cell culture freezing container from the -80°C freezer.
- Follow steps 1-2d above.
- Following trypsinization, resuspend cells in 9 mL of freezing medium containing 90% FBS and 10% DMSO.
- Aliquot 1 mL of re-suspended cells into each cryovial.
- Place cryovials containing cells into the pre-chilled cell culture freezer and replace into the -80°C freezer.
- Keep cells at -80°C overnight, and move to the liquid nitrogen tank the next day (tower 4, box 1).