

FOOD SAFETY LAB / MILK QUALITY IMPROVEMENT PROGRAM



Standard Operating Procedure

Title: Determination of Silage pH

SOP #: 7.25 Version: 02 Revision Date: 4-16-2020 Effective Date: 6/26/2020

Author: Nicole Martin, Rachel Evanowski Approved by:

Determination of Silage pH

FILE NAME: 7.25-Determination of Silage pH.docx









TABLE OF CONTENTS

1.	INTRODUCTION Purpose Scope Definitions Safety	3
2.	MATERIALS	4
3.	PROCEDURE	5
4.	TROUBLESHOOTING	6
5.	REFERENCES	6
6.	METHOD REVIEWS & CHANGES	6

Determination of Silage pH FSL/MQIP @ CORNELL UNIVERSITY

Effective 06/26/2020 Version 02 Revised 04/16/2020

SECTION 1 INTRODUCTION

1.1 Purpose

The purpose of this document is to set forth standard guidelines for performing pH measurements on silage.

1.2 Scope

This SOP applies to the Food Safety Lab and the Milk Quality Improvement Program. The protocols may also be used by laboratory members from other locations.

1.3 Definitions

DI water – Deionized water

Silage – Refers to all types of fermented animal feed including, but not limited to, corn silage and hay silage.

1.4 Safety

Wear gloves, safety glasses, and other appropriate personal protective equipment for the entire procedure.

Determination of Silage pH FSL/MQIP @ CORNELL UNIVERSITY 20 Version 02 Revised 04/16/2020

Effective 06/26/2020

SECTION 2 MATERIALS

- Balance
- DI Water
- pH Meter
- Stomacher and stomacher bags
- Beaker



Determination of Silage pHFSL/MQIP @ CORNELL UNIVERSITY

Effective 06/26/2020 Version 02 Revised 04/16/2020

SECTION 3 PROCEDURES

3.1 Sample Preparation

- 3.1.1. Weigh 5g of silage into a stomacher bag.
- 3.1.2. Add 45mL of DI water to the silage in the stomacher bag.
- 3.1.3. Stomach the silage and DI water at 230 rpm for 60 seconds.
- 3.1.4. Allow the stomached silage to sit undisturbed for 30 minutes prior to taking pH measurement

3.2 pH Measurement

3.2.1. Take pH measurement of silage slurry according to the pH meter manufacturer's instructions.



Determination of Silage pH FSL/MQIP @ CORNELL UNIVERSITY | O20 | Version | O2 | Revised | O4/16/2020

Effective 06/26/2020

SECTION 4 TROUBLESHOOTING

SECTION 5 REFERENCES

SECTION 6 METHOD VERSION & CHANGES

VERSION	DATE	EDITOR	COMMENTS
Version 1	11/5/2014	Nicole	Original SOP
Version 2	4/16/2020	Rachel E.	Formatting changes to fit new SOP format.
			Grammatical errors were fixed in the Procedures.
Version 3			