



FOOD SAFETY LAB / MILK QUALITY  
IMPROVEMENT PROGRAM

*Standard Operating Procedure*



Title: **actA PCR**

SOP #: 8.1.1.2.2.1

Version: **02**

Revision Date: **04/28/2020**

Effective Date: **05/11/2020**

Author: **Jordan Skeens, Emily Wright**

Approved by:

**actA PCR**

**FILE NAME: actAPCR.doc**



## TABLE OF CONTENTS

1.	INTRODUCTION	3
	Purpose	
	Scope	
	Definitions	
	Safety	
2.	MATERIALS	4
3.	PROCEDURE	5
4.	TROUBLESHOOTING	6
5.	REFERENCES	7
6.	METHOD REVIEWS & CHANGES	8



## SECTION 1 INTRODUCTION

### 1.1 Purpose

A standard *actA* hot start PCR method that allows differentiation of the *L. monocytogenes* control strain FSL R3-0001 from wild-type *L. monocytogenes* strains. FSL R3-0001 is currently used as a genetically modified control strain by Silliker laboratories and this procedure will allow the user to determine whether *L. monocytogenes* isolates from test samples represent wild-type isolates or contamination with FSL R3-0001.

This SOP is derived from an earlier protocol titled “*actA* PCR ASSAY TO DIFFERENTIATE LISTERIA MONOCYTOGENES FSL R3-001 CONTROL STRAIN FROM WILDTYPE STRAINS” authored by Emily Wright. The original protocol can be found archived on the Food Safety Wiki.

### 1.2 Scope

This SOP applies to the Food Safety and MQIP Labs.

### 1.3 Definitions

**BP:** base pair

**dNTPs:** deoxyribonucleotide monomers, which are the building blocks of DNA.

There are four dNTPs: cytosine, adenine, guanine and thymine

**GoTaq Green:** A PCR kit containing thermostable Taq polymerase, MgCl<sub>2</sub>, and a PCR buffer which contains a green loading dye. This loading dye does not interfere with the PCR and saves the step of loading dye incorporation prior to running gel electrophoresis.

**PCR:** Polymerase chain reaction is a process used to amplify a specific region of a DNA sequence.

**HOT START PCR:** A version of PCR (polymerase chain reaction) used to reduce nonspecific amplification and primer dimers.

### 1.4 Safety

PPE (gloves, lab coats and safety glasses) must be worn when working in the lab.



## SECTION 2 MATERIALS

### Single Tube Reactions

- Sterile water - Room 358B and Room 350A
- Sterile 0.2 mL tubes - Room 358B and Room 350A
- mpl-xbaF (forward) Primer, 12.5  $\mu$ M (5' - GCT CTA GAT CAA CAA GCA GCG AAA GATT - 3') - Room 352C, chest freezer
- actA MassR (reverse) Primer, 12.5  $\mu$ M (5' TTT ATG TGG TAA TTT GCT GTC 3') - Room 352C, chest freezer
- GoTaq Flexi Polymerase - Room 352, "Dumb" freezer
- GoTaq Green 5x PCR Buffer - Room 352, "Dumb" freezer
- dNTP solution consisting of 1 mM each dATP, dGTP, dCTP, dTTP. dNTP stocks in Room 352, "Dumber" freezer, in covered cryoblock
- MgCl<sub>2</sub>, 25 mM - Room 352, "Dumb" freezer
- Vortex Mixer - Room 358B
- Miniature Centrifuge - Room 358B
- Thermocycler - Room 356
- Micropipette and sterile filter tips
- Crushed Ice
- Controls: FSL R3-0001 and FSL X1-0001 (Both of these isolates must be used for positive controls to confirm amplicon size of unknown. FSL R3-0001 is the actA mutant, and FSL X1-0001 is the wild-type control, 10403S.

### High Throughput Reactions (96-well plate)

- Sterile reservoir basin - Room 358B
- Sterile 96-well plates - Room 358B
- Multi-channel pipette
- Sterile water - Room 358B and Room 350A
- mpl-xbaF (forward) Primer, 12.5  $\mu$ M (5' - GCT CTA GAT CAA CAA GCA GCG AAA GATT - 3') - Room 352C, chest freezer
- actA MassR (reverse) Primer, 12.5  $\mu$ M (5' TTT ATG TGG TAA TTT GCT GTC 3') - Room 352C, chest freezer
- GoTaq Flexi Polymerase - Room 352, "Dumb" freezer
- GoTaq Green 5x PCR Buffer - Room 352, "Dumb" freezer
- dNTP solution consisting of 1 mM each dATP, dGTP, dCTP, dTTP. dNTP stocks in Room 352, "Dumber" freezer, in covered cryoblock
- MgCl<sub>2</sub>, 25 mM - Room 352, "Dumb" freezer
- Vortex Mixer - Room 358B
- Miniature Centrifuge - Room 358B
- Thermocycler - Room 356
- Adhesive aluminum foil lids - Room 358B
- Micropipette and sterile filter tips
- Crushed Ice



## SECTION 3 PROCEDURES

### 3.1 Lysate Preparation

Refer to the “Preparing Cell Lysates for PCR” SOP on the Food Safety Lab wiki.

#### PCR Amplification

3.1.1 To prepare 200 $\mu$ L of the 1mM dNTP working solution:

3.1.1.1 Always make at least 200  $\mu$ L to avoid pipetting 1  $\mu$ L volumes, which can be inaccurate.

3.1.1.2 Thaw the individual 100mM dNTP tubes on ice.

3.1.1.3 Fill a 1.5ml tube with 192  $\mu$ L of dH<sub>2</sub>O.

3.1.1.4 Transfer 2  $\mu$ L of each dNTP into the 1.5mL tube and gently vortex to mix.

3.1.1.5 If only small portions of the 200 $\mu$ L will be used in your PCRs, aliquot this dNTP solution into 1.5ml tubes at appropriate volumes for your needs to reduce freeze thaw cycles.

3.1.1.6 This solution should be frozen at -20°C after use. 3.2.1.7

3.1.2 To prepare 100  $\mu$ L of the 12.5uM primer working solutions:

3.1.2.1 If the primer stock has not been reconstituted, refer to the “8.1.1.1.7-Primer Ordering and Reconstitution” SOP on the Food Safety Lab wiki.

3.1.2.2 Thaw the 100 $\mu$ M primer stock solutions on ice.

3.1.2.3 Fill two 1.5ml tubes with 87.5ul of 10mM Tris-HCL(pH8.0).

3.1.2.4 Add 12.5  $\mu$ L of your forward primer to one of the 1.5ml tubes, and 12.5  $\mu$ L of your reverse primer to the other 1.5ml tube. Label the tubes appropriately.

3.1.2.5 Gently vortex each tube to mix.

3.1.2.6 These solutions should be frozen at -20°C after use.

3.1.3 Working in Room 358B, prepare a master mix and Taq mix of the components listed in Table 1.

3.1.3.1 **\*\*IMPORTANT\*\***: Do not vortex the stock Taq polymerase solution or the Taq mix after the Taq is added, it will be damaged by vortexing.

3.1.3.2 Calculate master mix volumes required to include positive and



negative controls and 1-2 additional reactions to compensate for pipetting errors.

- 3.1.3.3 Store all master mix components and the prepared master mix on crushed ice.
- 3.1.3.4 Vortex and centrifuge down the contents of all reagent containers before opening (except Taq)
- 3.1.3.5 Pipette up and down to mix after each addition and to ensure that no reagent clings to the



3.1.3.6 pipette tip.

**Table 1: *actA* PCR Master Mix**

<u>Master Mix:</u>	<u>Volume per reaction</u>	<u>Final reaction concentration</u>
• dH <sub>2</sub> O	18.0 μL	
• 5X PCR buffer	8.0 μL	
• MgCl <sub>2</sub> (25mM)	3.0 μL	1.5mM
• dNTPs (1mM each)	5.0 μL	100uM
• mpl-xbaF (12.5μM)	2.0 μL	0.5uM
• actA MassR (12.5μM)	2.0 μL	0.5uM
	-----	
	38 μL per tube	
<u>Taq Mix:</u>	<u>Volume per reaction</u>	
• dH <sub>2</sub> O	7.6 μL	
• 5X PCR buffer	2.0 μL	
• Taq polymerase (GoTaq)	0.4 μL	
	-----	
	10 μL per tube	

3.1.4 Manual Hot Start

3.1.4.1 Briefly shake the tube containing the prepared master mix, but do not vortex

3.1.4.2 Spin briefly in the mini centrifuge to concentrate all contents at the bottom of the tube.

3.1.4.3 Dispense 38 μL of prepared master mix (but not the Taq Mix) into each well (one well per lysate to be amplified) in a 96-well PCR plate or into individual 0.2 mL PCR tubes. Cover the plate loosely or cover the tubes and hold on ice.

3.1.4.4 Remove the plate or tubes containing the master mix from room 358B to a lab bench.

3.1.4.5 Briefly spin down the previously-prepared lysates in a centrifuge, then add 2μL of lysate to its corresponding PCR well or tube. Pipette up and down to mix.

3.1.4.6 Briefly centrifuge the PCR plate or tubes to concentrate reagents at the bottom of the tube.



3.1.4.7 Perform the reaction using the following set of cycling conditions:

Hot Start PCR Conditions:

94°C	3 min	40 cycles
94°C	1 min	
53°C	1 min	
72°C	2 min	
72°C	5 min	
4°C	Hold	

3.1.4.8 Place the rack in the PCR machine, with a compression pad, and begin the cycle. Once the sample completes its initial (3 min at 94°C) and first denaturation step (1min at 94°C) and the temperature begins to decrease, pause the program when it reaches 80°C, and wait for the temperature to stabilize.

3.1.4.9 ***Be sure to thoroughly mix the Taq mixture before adding to the PCR tubes and always keep Taq on ice.*** Add 10µL of the Taq Mix to each PCR tube. Ensure all the tube caps are tightly closed, then close the thermocycler lid, wait 1 min, and resume the program.

3.1.5 Promptly remove the reactions from the thermocycler when complete. PCR product can be stored at 4°C for up to 1 week or at -20°C for up to 6 months.

**3.2 Gel Electrophoresis Confirmation of Amplification**

Refer to the “Gel Electrophoresis of PCR Product” SOP on the Food Safety Lab wiki.





### **3.3 Confirmation of results:**

3.4.1 FSL R3-0001 produces an amplicon of 600bp. FSL X1-0001 produces an amplicon of 2400bp. Compare the results of your unknown to these controls to determine whether it has the *actA* deletion, or it is a wild-type.

Refer to the “PCR Product Purification and Sanger Sequencing Submission” SOP on the Food Safety Lab wiki.



## SECTION 4                      TROUBLESHOOTING

### 4.1      **Issues with MgCl<sub>2</sub>**

- 4.1.1      MgCl<sub>2</sub> forms a concentration gradient when frozen and needs to be vortexed prior to use.
- 4.1.2      Every PCR reaction has an optimal MgCl<sub>2</sub> concentration range, usually between 1-4 mM. Since Mg<sup>+2</sup> forms complexes with dNTPs and acts as a cofactor for polymerase you may need to try several conditions to optimize your concentration.

### 4.2      **Too Much Enzyme**

- 4.2.1      Excess enzyme in your PCR can lead to smearing of PCR products.

### 4.3      **Wrong Primer Concentration**

- 4.3.1      If you have too little primer you won't see any product.
- 4.3.2      Too much primer may result in primer dimerization, which may reduce amplification.

### 4.4      **Wrong PCR Program**

- 4.4.1      Check your program while it's cycling to make sure it is the right program.

### 4.5      **Excess or Insufficient Template**

- 4.5.1      Too much template can inhibit PCR by binding all of the primers. If the lysate appears very turbid, dilute the template 1:10 or 1:100 and then use 2 μL of the diluted template for the PCR.
- 4.5.2      Too little template and amplification may be undetectable. Make sure you correctly prepared your dirty lysate and that you took a sufficiently large sample from your culture.

### 4.6      **Too Much dNTP or Degraded dNTP**

- 4.6.1      Excess dNTP inhibits PCR due to MgCl<sub>2</sub> depletion. Between 40-200 μM final concentration is optimal. dNTPs are sensitive to repeated freeze-thaw cycles. Make small aliquots when you get a fresh batch and turn over your stock frequently.
- 4.6.2      Note that the dNTP solution concentration suggested by the GoTaq Flexi manufacturer Promega suggests having each dNTP at 10mM, while this protocol uses 1mM due to less MgCl<sub>2</sub> being used.
- 4.6.3      The 1mM dNTP concentration has been validated for this PCR protocol, but should not be transferred for use in other PCRs, especially if the fragment length is much longer



## SECTION 5 REFERENCES

Frame, Peter. (2010). Ten Things That Can Kill Your PCR. *Bio-Synthesis*. [www.biosyn.com](http://www.biosyn.com).

## SECTION 6 METHOD VERSION & CHANGES

VERSION	DATE	EDITOR	COMMENTS
Version 1	01/30/2020		Original SOP
Version 2	04/28/2020	Ser15	Reformatted to new SOP template