

Gram-negative Genomic DNA Extraction (Qiagen):

Materials	QiaAMP DNA Mini kit: <ul style="list-style-type: none"> • Buffer ATL • Proteinase K solution • RNase • AL Buffer • AW1 Buffer • AW2 Buffer 10 mM Tris-HCl 100% Ethanol 2 heat blocks
Notes	Step
	Place 1.8ml of bacterial suspension into 2 ml tubes
	Find all heat blocks necessary and pre-heat. ____ 56C heat block ____ 70C heat block
	Prepare 1.5mL of 10 mM Tris-HCl in a 2 ml tube: ____ 15 ul 1 M Tris-HCl ____ 1485 ul dH ₂ O
	____ Pellet bacteria in microcentrifuge tube ____ Centrifuge 13,000 rpm for 7 minutes
	____ Remove supernatant from pellet
	____ Add 180 ul of Buffer ATL to the tubes ____ Pipet to homogenize the suspension
	____ Add 20 ul proteinase K and mix by pulse-vortexing for 15 s.
	Incubate at 56C for 45-60 minutes (depending on organism); flick tubes or vortex (to mix solution) every 15 minutes of incubation Time in: ____ Time: ____ Flick# ____ (Time: ____ Flick# ____) Time: ____ Flick# ____ Time out: ____ ____ Centrifuge briefly to collect liquid
	____ Add 4 ul RNase (100 mg/ml) ____ Gently mix by pulse-vortexing for 15 s ____ Centrifuge for a few seconds to collect liquid Incubate at room temperature for 10 minutes Start incubation: ____ End incubation: ____
	Pre-warm 10 mM Tris-HCL at 56C
	____ Add 200 ul Buffer AL , pulse-vortex for 15 s ____ Centrifuge for a few seconds to collect liquid

	<p>Incubate at 70C for 10 minutes Time in heating block: _____ Time out of heating block: _____</p>
	<p>____ Add 200 ul ethanol (96-100%) and pulse-vortex for 15 s. ____ Centrifuge for a few seconds to collect drops from inside the lid</p>
	<p>____ Add mixture to spin column in a 2 ml collection tube ____ Centrifuge at full speed (20,000 x g) for 1 minute. When centrifuging each wash, switch the direction of the caps. Tab direction?: _____</p>
	<p>____ Place the mini spin column in a clean 2 ml collection tube Discard the tube containing the filtrate.</p>
	<p>____ Open the spin column. Add 500 ul Buffer AW1. Close the cap. ____ Centrifuge at full speed for 1 minute. Tab direction?: _____</p>
	<p>____ Place the mini spin column in a clean 2 ml collection tube Discard the tube containing the filtrate.</p>
	<p>____ Open the spin column. Add 500 ul Buffer AW2. Close the cap. ____ Centrifuge at full speed for 3 minutes. Tab direction?: _____</p>
	<p>____ Place the mini spin column in a clean 2 ml collection tube Discard the tube containing the filtrate.</p>
	<p>____ Open the spin column. Add 500 ul Buffer AW2. Close the cap. ____ Centrifuge at full speed for 3 minutes. Tab direction?: _____</p>
	<p>____ Place the mini spin column in a clean 2 ml collection tube Discard the tube containing the filtrate.</p>
	<p>____ Run an empty spin at full speed for 1 minute. ____ Place the mini spin column in a clean 2 ml collection tube Discard the tube containing the filtrate.</p>
	<p>____ Run an empty spin at full speed for 2 minutes. ____ Place the mini spin column in a clean 1.5 ml collection tube (labeled elution1) Discard the tube containing the filtrate.</p>
	<p>____ Add 50 ul 10 mM Tris-HCl directly to the column. ____ Incubate at room temperature for 4 minutes Time start: _____ Time stop: _____</p>
	<p>____ Centrifuge at 8000 rpm for 1 min</p>
	<p>____ Place the mini spin column in a clean 1.5 ml collection tube (labeled elution2)</p>
	<p>____ Add 50 ul 10 mM Tris-HCl directly to the column. ____ Incubate at room temperature for 4 minutes Time start: _____ Time stop: _____</p>
	<p>____ Centrifuge at 8000 rpm for 1 min</p>
	<p>____ Nanodrop, then store both elutions at -20C</p>