

Tissue culture assay to quantify bacterial adhesion to HeLa cells

Created: July 18, 2000

Created by: Mario Roma and Martin Wiedmann

Last Revision: October 15, 2012

- 1) Grow HeLa cells to confluence in T-75 (approximately 8.0×10^6 cells)
- 2) Remove media, wash cell monolayer 3X with PBS
- 3) Remove PBS, add 1ml Trypsin-EDTA
- 4) After cells have detached, add 2ml DMEM+ 10% FBS w/o antibiotics (This should give approximately 2.67×10^6 cells/ml).
- 5) Pipette 9ml DMEM+ 10%FBS, w/o antibiotics into a 15ml centrifuge tube.
- 6) Add 1ml of cell suspension. (This should yield a total of 10ml at 2.67×10^5 cell/ml).
- 7) Plate 5 wells of a 6 well plate (with 4 sterile 12mm glass cover slips in each well) at 2ml per well. At this cell density, the plate should be ready to use the next day.
- 8) Pipette 19ml DMEM+ 10%FBS, w/o antibiotics into a 50ml centrifuge tube.
- 9) Add 1ml of cell suspension from step 4. (This should yield a total of 20ml at 1.3×10^5 cells/ml).
- 10) Plate 20 wells of a 24 well plate (with 1 sterile 12mm glass cover slip in each well) with 1ml per well. At this cell density, the plate should be ready to use the next day.
- 11) Grow overnight cultures of *E. coli* and *Erwinia* (or other bacterial species of interest).
- 12) Pipette 250 μ l of turbid broth (overnight culture) into 1.5ml centrifuge tubes.
- 13) Centrifuge at 13,000rpm for 5 minutes.
- 14) Remove supernatant, and resuspend pellet in 250 μ l of PBS+Mg.
- 15) Inoculate each well of both the 6 and 24 well plates with 50 μ l (for 6 well plate) and 10 μ l (for 24 well plate), respectively, of cell suspension.
- 16) Serially dilute cell suspension in PBS to 10^{-6} , and plate 100 μ l of the 10^{-5} and 10^{-6} dilutions on BHI agar plates. (i.e., a 10^{-6} and 10^{-7} final dilution).
- 17) Put the BHI agar plates in a 37 $^{\circ}$ C incubator overnight.
- 18) Place the 6 and 24 well plates in a 37 $^{\circ}$ C incubator with CO $_2$ for two hours.
- 19) After 2 hours, remove plates from incubator, remove media from each well, and wash each well containing cover slips 3X with PBS+Mg. After the third wash, add PBS+Mg to each well.

20) Remove cover slips and place in 15ml centrifuge tubes containing 5ml of .1% Triton X100. For the 6 well plates, all 4 coverslips from each well in the 6 well plate are added to one 15 ml tube. For the 24 well plate, all 4 coverslips from the 4 wells inoculated with the same bacterial strains are added to a single 15 ml tube. Thus, each 15 ml tube contains 4 coverslips.

21) Vortex each tube for 30 seconds.

22) Serially dilute the cell lysate to 10^{-3} , then plate 100 μ l of each dilution on BHI agar plates. Plate the 10^0 dilution once, plate the other dilutions in duplicate.

23) Put plates in a 37 $^{\circ}$ C incubator overnight.

24) Count colonies and tabulate results.