Evolution of natriuretic peptides-Current applications in human and animal medicine

The ABC of natriuretic peptides

- Atrial natriuretic peptide (ANP)
 - Produced in both atria
 - Triggered by increased atrial-wall tension following increased intravascular volume or pressure
 - Stored in granules to be released
 - Urodilatin-version of ANP found in the kidneys
 - Produced in renal distal tubules
 - Locally controls sodium and volume
 - Decreases renin secretion from macula densa
 - Inhibits aldosterone release from zona glomeruosa
 - Inhibits angiotensin II-stimulated proximal tubular Na⁺ and H₂O reabsorption
- B-type natriuretic peptide (BNP)
 - Found predominantly in the ventricles
 - Rapidly produced by cardiomyocytes w/ myocardial stretch & hypoxia
 - Actions:

| • | Diuresis and natriuresis via direct tubular |
|---|---|
| | actions (as above for urodilatin) |

- Inhibits cardiac sympathetic nervous system
- Inhibits RAAS by suppressing norepinephrine and aldosterone levels
- Anti-proliferative, antifibrotic and relaxation effects on heart and vascular tissues
- C-type natriuretic peptide (CNP)
 - Paracrine and autocrine functions
 - Found in vascular endothelial cells, central nervous system, kidney, chondrocytes, pituitary gland

Renal

↑GFR

♦Na+ resorbtion

Vascular

♦Arterial Tone


↓Venous tone

Antiproliferative

- Actions:
 - Neurotransmitter function
 - Regulation of vascular tone
 - Inhibition of fibroblast activation
 - · Role in enchochondral ossification
 - Role in testicular endocrine function, spermatogenesis and adequate placental function during pregnancy
 - Minimal diuretic/natriuretic actions
- Dendroaspis natriuretic peptide (DNP)
 - Isolated from the venom gland of the green mamba-unknown function
- Ventricular natriuretic peptide (VNP)
 - Role in the fluid homeostasis of primitive ray-finned bony fish

Natriuretic peptide receptors and neutral endopeptidase

- Natriuretic peptide receptor A
 - Expressed in lungs, brain, heart, adrenals, kidney, and terminal ileum
 - Actions:
 - Production of intracellular cGMP
 - Leads to further signal transduction
 - Natriuresis, inhibition of renin and aldosterone
 - Vasorelaxation
 - Anti-fibrotic, anti-hypertrophic and lusitropic
- Natriuretic peptide receptor B
 - Expressed in lungs, brain, skin, adrenals, kidneys, uterus and ovaries
 - Expression in veins predominantly over arteries
 - Activated predominantly by CNP
 - Actions:
 - Production of intracellular cGMP
 - Mediates vasorelaxation by actions on smooth muscle
 - Has a hypotensive effect in comparison to NPR-A



Cardiac

Lusitropic

Antifibrotic

Antiproliferative




SNS/RAAS

↑Vagal tone

♥SNS activity

♥Renin release

♦Aldosterone release



- Natriuretic peptide receptor C
 - · Expressed in lungs, brain, heart, adrenals, kidneys, mesenterium, fat tissue, placenta, veins and aorta
 - Modulates naturetic peptide activity at target organs-ie downregulates NP activity
- Natriuretic peptide receptor D-only found in the eel
- Neutral endopeptidase
 - Induced by adenylate cyclase, glucocorticoid, thrombin, calcitonin, and cytokines
 - Actions
 - Inactivation of enkephalins, neuropeptides (ie substance P) and β -amyloid peptide
 - Active elimination of natriuretic peptides from circulation by hydrolysis

Processing of pre-proBNP and its products

- Pre-proBNP
 - Stimulated production by cardiomyocytal stretch and ischemia/hypoxia
 - Also stimulated by
 - Endothelin-1
 - Angiotensin-II
 - Interleukin 1 β
 - Adrenergic agonists
 - Rapid removal of an amino acid terminal results in proBNP
- proBNP cleaved into two parts by corin or furin
 - Forms NT-proBNP and the active substance BNP
- Elimination of BNP via NPR-C, neutral endopeptidase activity, and renal clearance
- NT-proBNP cleared entirely by the kidneys

BNP and NT-proBNP testing in humans and animals

- Human testing considerations
 - BNP levels are decreased in obesity and increased in kidney disease
 - Excellent in predicting prognosis in short and long term for acute or chronic heart failure
 - Also useful in the full spectrum of cardiovascular diseases
 - Decreases in NT-proBNP was associated with improved outcome
- Veterinary applications
 - BNP able to discriminate between heart failure and non-heart failure in dogs and is correlated with severity of heart failure
 - Confirmed on two separate studies (DeFrancesco 2007 and Oyama 2008)
 - In comparing BNP to NT-proBNP, the later was diagnostically superior in evaluating canine respiratory distress
 - Aortic banding in 6 male beagles showed NT-proBNP to be better correlated with LVEDP compared to NT-proANP

Ouestions

- 1. Interaction of BNP with NPR-A will ultimately result in
 - a. Inositol triphosphate production
 - b. Protein kinase activation
 - c. Diacylglycerol production
 - d. Adenylyl cyclase activation
- List two mechanisms of elimination of BNP
 - a. NPR-C
 - b. Neutral endopeptidases
- 3. As opposed to ANP, BNP is
 - a. Stored in granules for release
 - b. Triggered by increased atrial wall tension
 - c. Homologous to urodilatin
 - d. Stimulated for production during stretch of cardiomyocytes
- 4. CNP predominantly exerts its effects
 - a. In a paracrine and autocrine manner
 - b. On angiotensin II-mediated sodium reabsorption mechanisms
 - c. On aquaporin channels in ray-finned bony fish
 - d. On osteoblasts
- Diagram the processing of pre-proBNP to active BNP



