

2019-05-15 Meeting notes

Date

15 May 2019

Attendees

- Clayton L Birkett

Goals

Discussion items

Time	Item	Who	Notes
	GrainGenes update	Taner	
	T3	Clay	1000 exomes genotype experiment page select lines link to download vcf
	T3	Vic	1000 exomes passport data
	T3	Clay	primer selection shiny-server http://tcapg.ag.cornell.edu/primer_filter/
	BreedBase		2011-13 genotype experiments loaded
	Wheat PHG	Clay	options for reference intervals • what is minimum size? 25 bases • which intervals to combine? 10,000 bases GBS + Exome + gene model - 250,000 reference ranges Exome + gene model - 120,000 reference ranges TASSEL functions exclude gene models if no exome or GBS sequence
	GOBii	Clay	progress on new install

Action items

- ☒ Clayton L Birkett- ask Lukas if wizard can work on projects